Maximum convenience and great fun under the shower and on the toilet with the Dukki and the Turtel







Huka



For many years already, Huka have been specialists in the field of children's rehabilitation. We invent and produce our products in the Netherlands, together with designers and rehab specialists. New products are often created in collaboration with children because they know best what's good for them. Our products are tested extensively because, in terms of quality, only the best of the best is good enough for your children. Our goal: we want children, their parents, and their carers to consider our products to be welcome and excellent support aids for sitting and moving around.



Children with a handicap often need help showering or going to the toilet. In this brochure it gives us great pleasure to present our shower and commode chairs, Turtel and Dukki, products that make these daily necessities easier and more fun. Sitting comfort, dimensions, hygiene,

ease of use, and design all played an important role in the development of these two chairs.

Not only your child will benefit. You will also find out how convenient our products will be to you. You will have both hands free to help your child on

the toilet or in the shower. There is no need for heavy lifting to support your child, and when your child is finished, you simply toll the chair aside. On the next pages you can see just how wonderful and useful our Turtel and Dukki are.





Benefit from the convenience of our Turtel on the toilet





Accessories

Several accessories are available. A footrest will enable you child to rest its feet. An inlay cushion will provide an even more comfortable position to sit. Of course padded lumbar, chest and feet straps can be added for maximum support. A headrest can be fitted to support the head. You can also order a collection pan for underneath the seat.

Easy to use and comfortable

Standard features of the Turtel include a plastic bucket-shaped seat, a galvanized steel frame and plastic castor wheels. The seat shape guarantees perfect support of your child's upper body. The seat height can be adjusted easily by adjusting the legs. Its reach is such that the Turtel also fits perfectly on top of a raised toilet. The adjustable legs also allow you to tilt back the Turtel into several positions. Due to the unique

wheels will continue to provide excellent manoeuvrability in any tilting position. For safety reasons, the castors have been fitted with brakes in the front. This prevents the chair from rolling away accidentally

mounting system of the castors, the

Smart dimensions

during use.

The chair is available in no less than three sizes, which makes it suitable for children from the ages of 4 to 12.

For older children, the Turtel Supermaxi is available.

Depending on its size, the Turtel is available with a raspberry red, capri blue, or water green bucket seat and black frame. All sizes are also available with a neutral light-grey seat with white frame.





Double convenience with our Dukki, the 2-in-1 chair for use in the toilet and the shower

The Dukki is a multi-functional toilet/shower chair. The chair is easy to use and easy to maintain. The material is cleaned with a simple washing.

Smart dimensions

The Dukki is available in no less than four sizes for children from the ages of 4 to 15. Since every child is different, all Dukki sizes can be fitted with two different sizes abduction blocks and armrests. The selection of a cushion allows further customisation of the dimensions to your child. This is how we achieve maximum comfort and good support.

All sizes are available in two colours: light blue and light green.

Very comfortable

Standard features of the Dukki include an aluminium frame with plastic castor wheels, a plastic bucket-shaped seat, and a plastic footrest. The seat has been designed to ensure that your child is sitting comfortably. The seat can be tilted easily for use in the shower. Your child can relax and fully enjoy the warm water.

Accessories

Several accessories are available including a padded cushion, an abduction block with armrest. padded lumbar, chest and feet straps, and a headrest. You can also order a collection pan for underneath the seat.

Great ease of use

The Dukki will serve you many times a day. You can use it fast and with ease. The armrest and the abduction block can be removed in one simple movement. This will enable you to lift your child into and out of the chair easily. In addition to its regular function to rest the arms, you can also use the armrest to put a toy on it to divert your child's attention during washing.

The tilt position is ideal to position your child correctly in the chair when using a lift, for instance. This will often be the case when your child is growing and it becomes too much of a strain to lift the child up by hand. Both the tilt angle of the seat and the headrest height can be adjusted without any tools. The same applies

to the footrest.

The supplied tool can be used to adjust the frame height, e.g. to the height of your toilet.

You can also use this tool to disassemble the Dukki fast and without any effort. This will enable you to take the Dukki along so it can serve you on your holiday as well.









specia		